



Letter from the President

Because of the unusually mild winter on Cape Cod this year, it seems like our laboratory should already be full of horseshoes crabs. Fortunately, winter allows us some time to squeeze in workshops or attend to writing and research.

Recently I discovered that since our founding, ACC has authored or co-authored over 100 original articles on LAL or endotoxin. And this does not include the numerous articles which have appeared in the LAL UPDATE over the past 15 years! Although ACC does not manufacture therapeutic drugs, devices (other than LAL), or reagents for tissue culture, we have tested these products or used them in our research. We therefore write from experience. In this LAL UPDATE, Dr. Dawson continues to provide practical information with an article that should be of special interest to cell culture and biotechnology specialists.

Remember our technical staff is always ready to help with any endotoxin testing/removal applications you may have. Please contact us by phone, fax, or e-mail and be sure to visit our web site.

Sincerely,

Thomas J. Novitsky, Ph.D.

The Significance of Endotoxin to Cell Culture and Biotechnology

By Michael E. Dawson, Ph.D.

Manufacturers of biopharmaceuticals have two concerns about endotoxin. First, the concentration of endotoxin in finished injectable products must be below the endotoxin limit. This concern applies to all injectable drugs and biologicals and non-pyrogenic medical devices. It is well documented and will not be discussed in this article. The second concern is the influences of endotoxin on the expression systems used to produce biopharmaceuticals and upon the products themselves.

The effects of endotoxin vary greatly for different cell types or cell lines. One should be aware of the effect upon the cells of interest so the appropriate level of control can be instituted. Even parent and daughter cell lines may show very different sensitivities to endotoxin.¹ Awareness of the surface receptors of the cells being used can alert an investigator to the potential for endotoxin sensitivity. The presence of the receptor, CD14, on human cells is strongly associated with endotoxin sensitivity.² This is acutely important where the culture medium contains serum. Serum components, such as lipopolysaccharide binding protein (LBP) and septicin complex, can potentiate endotoxin activation of CD14-bearing cells.^{3,4} Conversely, some mammalian and invertebrate cell systems are tolerant of endotoxin. The presence of endotoxin in the media for bacterial and yeast culture is generally less of an issue. The effects of endotoxin on examples of these biological systems considered below.

In studies involving cell culture, one must guard against spurious results caused by biological responses to contaminants. Also, the properties of the products must be distinguished from those of contaminant endotoxin. This is especially true for cytokines.

Membrane and Morphological Effects

Endotoxin is an amphipathic membrane component of gram negative bacteria. It is not surprising that it might interfere with the membrane structure and function in other cells. Through its phosphate groups it reacts with cationic proteins, and its fatty acid chains interact with lipid membranes and hydrophobic regions of proteins. The interaction of endotoxin with cell membranes may be apparent as morphological changes, such as perturbed Chinese hamster ovary cell membranes,⁵ surface ruffles and increased organelles,⁶ large vacuoles in the cytoplasm,⁷ morphological damage and decreased hexose uptake,⁸ and severe membrane damage.⁹ Sometimes morphological changes are reversible upon exposure to clean media.¹⁰ Also, the effects upon morphology may vary with different endotoxins.¹¹

If severe, the effects of endotoxin upon membrane structure and function may be cytotoxic. The degree of toxicity varies between different endotoxins¹² and between cell lines.¹³ With at least some cells, however, it is possible to induce tolerance to endotoxin.¹⁴

Secretion/Production

Membranes are critical to the synthesis of cell products and endotoxin can influence the processes. This is clearly of critical concern in cell based research and in the development of commercial cell lines. The sensitivities of cells are diverse; in one example, two non-responsive T cell lines were cloned from an endotoxin responsive parent.¹⁵ Also, the effects of endotoxin may be modified by other factors, such as temperature.¹⁶

The influence of endotoxin upon cytokine production by single cell types *in vitro* has been widely reported. Endotoxin has been cited as the most potent stimulus for tumor necrosis factor (TNF) production.¹⁷ Effects upon transcription have been noted^{18,19,20,21} and it has been suggested that the primary regulation of stimulation may occur at this level.²²

A range of cell products other than cytokines are also influenced by endotoxin. Products for which synthesis by mammalian cells is stimulated or enhanced include prostaglandin,²³ acid phosphatase,⁶ fibrinolytic inhibitor,²⁴ collagenase production,²⁵ nerve growth factor,²⁶ secretion of factor B,²⁷ polypeptides,²⁸ platelet activating factor,²⁹ adhesion inhibitor,³⁰ adhesion molecule-1³¹ and procoagulant (dose dependent).³²

Inhibitory effects have been observed for angiotensin converting enzyme activity,³³ and synthesis of proteoglycan²³ and alpha2 macroglobulin.³⁴ These responses may be dose dependent and in some cases are quite different for different endotoxins. For example, at some doses *Bacteroides* endotoxin stimulated phosphatase production, while that from *E. coli* was inhibitory.¹¹

Mitogenicity

Endotoxins are powerful mitogens to some mammalian cell types and lines. This effect can create problems in the culture of susceptible cells. On the other hand, some long-established cell lines are tolerant to high endotoxin concentrations in the culture medium. A correlation between mitogenicity and both LAL reactivity and pyrogenicity of several different endotoxins has been reported.³⁵ However, positive correlations are not always evident and endotoxins of equal mitogenicity may differ significantly in their LAL reactivity and pyrogenicity.³⁶ Therefore, LAL activity is not a good predictor of endotoxin mitogenicity.³⁷

Some mitogens are heavily contaminated by endotoxin (Table 1) and the high concentrations raise questions about the contribution of endotoxin to their mitogenicity. Similarities in the effects of endotoxin and mitogens have been reported.³⁸

Concanavalin A	<0.003 – 9.2 x 10 ⁵ ng/mg
Pokeweed mitogen	1.1 x 10 ⁶ ng/mg protein (<i>sic.</i>)

Table 1. Endotoxin concentrations in two mitogens³⁹

To add to the complexity, endotoxin can inhibit cell division¹⁰ and, as with stimulation of mitosis, these effects depend on the origin of the endotoxin.⁴⁰ There may be a concentration dependent effect, ranging from promotion of mitosis to inhibition at higher concentrations. It is interesting to note that the concentrations of endotoxin in some commonly used mitogens are as high as those found to inhibit cell division in some cell lines.

Other Effects

Other *in vitro* cellular effects of endotoxins reported in the literature include: influence upon adherence,⁴¹ decreased insulin binding and endocytosis,⁴² loss of TNF binding sites,⁴³ protection from HIV

infection,⁴⁴ tumoricidal activity,⁴⁵ and increases in potassium current and in the number of potassium channels.⁴⁶

Multiple Effects and Synergism

The multiple biological responses elicited by endotoxin may result in confusion about cause and effect. The influence of endotoxin upon secretion of cell products can be indirect. For example, an antibody to interferon was shown to prevent the apparent inhibition of elastase production by endotoxin. This suggests that the interferon produced by the cell caused the inhibition. In the presence of the antibody, endotoxin actually increased elastase secretion. The inhibitory effect was therefore secondary and due to endotoxin stimulated interferon production.⁴⁷ Similarly, an 80% inhibition of lipase production by cultured heart cells is mediated by endotoxin induced TNF produced by the cells.⁴⁸ These observations are a clear warning that observed effects may be indirect consequences of endotoxin exposure. Similarly, mitogenicity may be a secondary effect. Endotoxin stimulated IL-1 production by macrophages and neutrophils has been shown to be mitogenic in spleen cells.⁴⁹ In this case the effect was elicited by *Bacteroides gingivalis* endotoxin but not by that from enteric bacteria.

An interesting synergistic interaction between lymphokines and endotoxin results in the induction of chemiluminescence in macrophages. This response has been suggested as a replacement for a cytotoxicity test.⁵⁰ In another study, both synergism and antagonism between endotoxin and cytokines in the same system were reported.⁵¹ Granulocyte-macrophage colony stimulating factor makes cells responsive to endotoxin, which induces TNF production. TNF stimulates IL-6 production which inhibits further TNF production in a negative feedback loop.

In vitro Fertilization

In vitro fertilization is a special application of tissue culture technology. Improvements in the fertilization rate have been achieved in media of low endotoxin concentration.⁵² In the absence of detectable endotoxin in the media (at a detection limit of 0.1 ng/ml), a pregnancy rate of 27.2% was achieved, compared with only 6.3% when endotoxin was present. Morphological and ultrastructural changes were evident in the presence of endotoxin. In another study,⁵³ an increasing degree of fragments in the conceptus was observed at endotoxin concentrations above 1 ng/ml. The oocyte fertilization rate was 66% when endotoxin concentrations were less than 1 ng/ml in the medium, compared with 53% when greater than this. Differences were more marked when pregnancies were compared. The pregnancy rate was 8% if the endotoxin concentration in the medium exceeded 1 ng/ml but 32% when no endotoxin was detected (<0.1 ng/ml).

Conclusion

Because of the great variability in the response of cells to endotoxin, it is not possible to state a critical level at which endotoxin begins to interfere with cell function and growth, or what form that interference might take. The effects of endotoxin, and the sensitivity of the cells in question, can only be established by growing cells in a medium free of detectable endotoxin. This may or may not be practicable or possible. If it is, critical levels of contamination can be established and limits can be set.

The effects of endotoxin upon cellular processes has long been recognized. Particular attention has been paid to endotoxin contamination of sera⁵⁴ because of the widespread use of sera in cell culture media. Fumarola and Jirillo^{55,56} have called vigorously for recognition of the effect of endotoxins on cell lines and the need for appropriately low endotoxin concentrations. They state that "contamination with endotoxin of some preparations commonly used in biology may give false results." There has been at least one call for the "absolute absence" of endotoxin in tissue culture media.⁵⁷ These authors also discuss the use of endotoxin to stimulate cells and enhance their reactivity, obviously in a controlled way.

Endotoxin contamination of products should be minimized, and not only because of concerns about pyrogenicity and endotoxin limits. The properties of some products of biotechnology are shared by endotoxin. These products cannot be properly evaluated if the preparation is contaminated with endotoxins.⁵⁸ This was well illustrated by the experience of an investigator using an End-X B15 endotoxin affinity resin (Associates of Cape Cod, Inc.). Endotoxin was successfully removed from a protein solution to well below levels of concern. The specific enzymatic activity of the protein was unaffected by treatment. However, the mitogenicity of the preparation dropped dramatically. The investigator eventually concluded that the mitogenicity of the original preparation was probably due to the substantial endotoxin contamination and was not an inherent property of the protein.

The question of endotoxin contamination should be addressed early in product development and certainly before any scale-up. This will enable both minimization of deleterious effects and proper determination of the product's properties. Early attention to endotoxin related issues is then likely to follow through to full-scale production. This will have benefits in final release testing when the endotoxin specification is as important for biopharmaceuticals as it is for conventionally produced therapeutics.

References

- Ohki, K., A. Nagayama, and O. Kohashi. 1992. Differential effects of bacterial lipopolysaccharide and interferon-gamma on proliferation of two factor-dependent macrophage cell lines. *Cell Structure and Function* 17:161-167.
- Pugin, J., Schurer-Maly, C., Leturq, D., Moriarity, A., Ulevitch, R.J., and Tobias, P.S. 1993. Lipopolysaccharide activation of human endothelial and epithelial cells is mediated by lipopolysaccharide-binding protein and soluble CD14. *PNAS* 90:2744-2748.
- Mathison, J.C., Tobias, P.S., Wolfson, E., and Ulevitch, R.J. 1992. Plasma lipopolysaccharide (LPS)-binding protein. A key component in macrophage recognition of Gram-negative LPS. *J. Immunol.* 149:200-206.
- Wright, S.D., Ramos, R.A., Patel, M., and Miller, D.S. 1992. Septin: A factor in plasma that opsonizes lipopolysaccharide-bearing particles for recognition by CD14 on macrophages. *J. Exp. Med.* 176:719-727.
- Csaba, G., O. Torok, and P. Kovacs. 1984. Influence of endotoxin treatment on TSH binding and TSH-induced imprinting in CHO cell culture. *Cell Biol Int. Rep.* 8:1.
- Barker, K.A. and S. C. Holt. 1983. The *in vitro* effect of *Actinobacillus actinomycetemcomitans* strain Y4 lipopolysaccharide on murine peritoneal macrophages. *Can. J. Microbiol.* 29:1552-1563.
- Tagikuchi, S. and A. Koga. 1988. Effects of bile acids and endotoxin on the function and morphology of cultured hamster Kupffer cells. *Virchows Archives. B. Cell Pathology* 54:303-311.
- Pagani, R., M. T. Portolés, I. Díaz-Laviada, and A. M. Muncio. 1988. Morphological damage induced by *Escherichia coli* lipopolysaccharide in cultured hepatocytes: Localization and binding properties. *Br. J. Exp. Pathol.* 69:537-549.
- Paulsen, D.B., Mosier, D.A., Clinkenbeard, K.D., and Confer, A.W. 1989. Direct effects of *Pasteurella haemolytica* lipopolysaccharide on bovine pulmonary endothelial cells *in vitro*. *Am. J. Vet. Res.* 50:1633-1637.
- Fumarola, D., E. Jirillo, G. Miragliotta, and E. Magliulo. 1982. Influence of lipopolysaccharide from *Pasteurella multocida* on fibroblasts cultured *in vitro*. *IRCS. Med Sci.* 10:647-648.
- Johne, B., G. Gaudernack, and B. Morland. 1987. Effects of endotoxins from *Bacteroides intermedium* and *Escherichia coli* on human monocytes *in vitro*. *Acta Pathol. Microbiol. Immunol. Scand.* 95 (Suppl. C):241-250.
- Reitmeyer, J. C. and J. W. Peterson. 1988. *In vitro* cytotoxicity of lipopolysaccharides for chinese hamster ovary cells. *FEMS Microbiol. Lett.* 49:393-396.
- Harlan, J. M., L. A. Harker, G. E. Striker, and L. J. Weaver. 1983. Effects of lipopolysaccharide on human endothelial cells in culture. *Thromb. Res.* 29:15-26.
- Hartung, T. and A. Wendel. 1992. Endotoxin-inducible cytotoxicity in liver cell cultures—II. Demonstration of endotoxin-tolerance. *Biochem. Pharmacol.* 43:191-196.
- Mukaida, N., T. Kasahara, K. Shioiri-Nakano, and T. Kawai. 1984. Interleukin 2-dependent T-cell line acquires responsiveness to phorbol myristate acetate and lipopolysaccharide in the course of long term culture. *Immunol. Commun.* 13:475-486.
- Velasco, S., M. Tarlow, K. Olsen, J. W. Shay, G. H. McCracken, Jr., and P. D. Nisen. 1991. Temperature-dependent modulation of lipopolysaccharide-induced interleukin-1 and tumor necrosis factor expression in cultured human astroglial cells by dexamethasone and indomethacin. *J. Clin. Invest.* 87:1674-1680.
- Bienhoff, S. E., G. K. Allen, and J. N. Berg. 1992. Release of tumor necrosis factor-alpha from bovine alveolar macrophages stimulated with bovine respiratory viruses and bacterial endotoxins. *Veterinary Immunology and Immunopathology* 30: 341-357.
- Gessani, S., F. Belardelli, A. Pecorelli, P. Puddu, and C. Baglioni. 1989. Bacterial lipopolysaccharide and gamma interferon induce transcription of beta interferon mRNA and interferon secretion in murine macrophages. *Journal of Virology* 63:2785-2789.
- Zuckerman, S. H., G. F. Evans, Y. M. Snyder, and W. D. Roeder. 1989. Endotoxin-macrophage interaction: Post-translational regulation of tumor necrosis factor expression. *J. Immunol.* 143:1223-1227.
- Tannenbaum, C. S., T. J. Koerner, M. M. Jansen, and T. A. Hamilton. 1988. Characterization of lipopolysaccharide-induced macrophage gene expression. *J. Immunol.* 140(10):3640-3645.
- Ziegler-Heitbrock, H. W. L., M. Blumenstein, E. Käfferlein, D. Kieper, I. Petersmann, S. Endres, W. A. Flegel, H. Northoff, G. Riethmüller, and J. G. Haas. 1992. *In vitro* desensitization to lipopolysaccharide suppresses tumour necrosis factor, interleukin-1 and interleukin-6 gene expression in a similar fashion. *Immunology* 75:264-268.
- Taffer, S. M., K. J. Singhel, J. E. Overholzer, and S. A. Shurtleff. 1989. Regulation of tumor necrosis factor expression in a macrophage-like cell line by lipopolysaccharide and cyclic AMP. *Cell. Immunol.* 120:291-300.
- Morales, T. I., L. M. Wahl, and V. C. Hascall. 1984. The effect of bacterial lipopolysaccharides on the biosynthesis and release of proteoglycans from calf articular cartilage cultures. *J Biol Chem* 259:6720-6729.
- Dubor, E., A. M. Dosne, and L. A. Chedid. 1986. Effect of polymyxin B and colimycin on induction of plasminogen anti-activator by lipopolysaccharide in human endothelial cell culture. *Infect. Immun.* 52:725-729.
- Cury, J. D., E. J. Campbell, C. J. Lazarus, R. J. Albin, and H. G. Welgus. 1988. Selective up-regulation of human alveolar macrophage collagenase production by lipopolysaccharide and comparison to collagenase production by fibroblasts. *J. Immunol.* 141(12):4306-4312.
- Mallat, M., R. Houlgatte, P. Brachet, and A. Prochiantz. 1989. Lipopolysaccharide-stimulated rat brain macrophages release NGF *in vitro*. *Dev. Biol.* 133:309-311.
- Miyama, A., Y. Kawamoto, H. Ichikawa, K. Okamoto, S. Hara, and T. Inoue. 1980. Complement proteins and macrophages. II. The secretion of factor B by lipopolysaccharide stimulated macrophages. *Microbiol. Immunol.* 24:1223-1232.
- Hamilton, T. A., M. M. Jansen, S. D. Somers, and D. O. Adams. 1986. Effects of bacterial lipopolysaccharide on protein synthesis in murine peritoneal macrophages: relationship to activation for macrophage tumoricidal function. *J. Cell Physiol.* 128:9-17.
- Wang, J., M. Kester, and M. J. Dunn. 1988. The effects of endotoxin on platelet-activating factor synthesis in cultured rat glomerular mesangial cells. *Biochim. Biophys. Acta* 969:217-224.
- Wheeler, M. E., F. W. Luscinskas, M. P. Bevilacqua, and M. A. Gimbrone, Jr. 1988. Cultured human endothelial cells stimulated with cytokines or endotoxin produce an inhibitor of leukocyte adhesion. *J. Clin. Invest.* 82:1211-1218.
- Wong, D. and K. Dorovini-Zis. 1992. Upregulation of intercellular adhesion molecule-1 (ICAM-1) expression in primary cultures of human brain microvessel endothelial cells by cytokines and lipopolysaccharide. *Journal of Neuroimmunology* 39:11-22.
- Delvos, U., B. Janssen, and G. Muller-Berghaus. 1987. Effect of lipopolysaccharides on cultured human endothelial cells. Relationship between tissue factor activity and prostacyclin release. *Blut* 55:101-108.
- Watanabe, K., G. Lam, R. S. Keresztes, and E. A. Jaffe. 1992. Lipopolysaccharides decrease angiotensin converting enzyme activity expressed by cultured human endothelial cells. *J. Cell Physiol.* 150:433-439.
- Bauer, J., U. Ganter, and W. Gerok. 1988. Endotoxin abolishes the induction of 2-macroglobulin synthesis in cultured human monocytes indicating inhibition of the terminal monocyte maturation into macrophages. *Adv. Exp. Med. Biol.* 240:425-432.
- Elin, R. J., A. L. Sandberg, and D. L. Rosenreich. 1976. Comparison of the pyrogenicity, Limulus activity, mitogenicity and complement reactivity of several bacterial endotoxins and related compounds. *J Immunol.* 117:1238-1242.
- Butler, T. and G. Moller. 1977. Mitogenic response of mouse spleen cells and gelation of Limulus lysate by lipopolysaccharide of *Yersinia pestis* and evidence for neutralization of the lipopolysaccharide by polymyxin B. *Infect. Immun.* 18:400-404.
- Alving, C. R. and E. C. Richardson. 1984. Mitogenic activities of lipid A and liposome-associated lipid A: effects of epitope density. *Rev. Infect. Dis.* 6:493-496.
- Letwin, B. W. and F. W. Quimby. 1987. Effects of concanavalin A, phytohemagglutinin, pokeweed mitogen, and lipopolysaccharide on the replication and immunoglobulin synthesis by canine peripheral blood lymphocytes *in vitro*. *Immunol. Lett.* 79-85.
- Rickles, E. R., J. Levin, J. A. Hardin, C. E. Barr, and M. E. J. Conrad. 1977. Tissue factor generation by human mononuclear cells: Effects of endotoxin and dissociation of tissue factor generation from mitogenic response. *J Lab. Clin Med* 89:792-803.
- Wille, J. J., J. Park, and A. Elgavish. 1992. Effects of growth factors, hormones, bacterial lipopolysaccharides, and lipoteichoic acids on the clonal growth of normal ureteral epithelial cells in serum-free culture. *J. Cell Physiol.* 150:52-58.
- Edelman, J., C. Cardozo, and M. Lesser. 1989. Lipopolysaccharide stimulates alveolar macrophage adherence *in vivo* and *in vitro*. *Agents Actions* 26:287-291.
- Portolés, M. T., R. Pagani, M. J. Ainaga, I. DíazLaviada, and A. M. Muncio. 1989. Lipopolysaccharide-induced insulin resistance in monolayers of cultured hepatocytes. *Br. J. Exp. Pathol.* 70:199-205.

43. Ding, A. H., E. Sanchez, S. Srimal, and C. F. Nathan. 1989. Macrophages rapidly internalize their tumor necrosis factor receptors in response to bacterial lipopolysaccharide. *J. Biol. Chem.* **264**:3924-3929.
44. Kornbluth, R. S., P. S. Oh, J. R. Munis, P. H. Cleveland, and D. D. Richman. 1989. Interferons and bacterial lipopolysaccharide protect macrophages from productive infection by human immunodeficiency virus in vitro. *J. Exp. Med.* **169**:1137-1151.
45. van Bossuyt, H., C. Desmaretz, B. Rombaut, and E. Wisse. 1989. Response of cultured rat Kupffer cells to lipopolysaccharide. *Archives of Toxicology* **62**:316-324.
46. Nelson, D. J., B. Jow, and E. Jow. 1992. Lipopolysaccharide induction of outward potassium current expression in human monocyte-derived macrophages: Lack of correlation with secretion. *J. Membr. Biol.* **125**:207-218.
47. Dodon, M. D. and S. N. Vogel. 1985. Analysis of effects of lipopolysaccharide and interferon on murine macrophages: modulation of elastase secretion in vitro. *Infect. Immun.* **49**:709-714.
48. Friedman, G., R. Gallily, T. Chajek-Shaul, O. Stein, E. Shiloni, J. Etienne, and Y. Stein. 1988. Lipoprotein lipase in heart cell cultures is suppressed by bacterial lipopolysaccharide: an effect mediated by production of tumor necrosis factor. *Biochim. Biophys. Acta* **960**:220-228.
49. Isogai, H., E. Isogai, N. Fujii, K. Oguma, W. Kagota, and K. Takano. 1988. Histological changes and some in vitro biological activities induced by lipopolysaccharide from *Bacteroides gingivalis*. *Zentralblatt für Bakteriologie Mikrobiologie und Hygiene, Series A* **269**:64-77.
50. Muto, S., H. Ogino, and H. Yuki. 1989. Synergistic effect of lymphokines and lipopolysaccharides on macrophage chemiluminescence. *J. Immunol. Methods* **111**:51-57.
51. Aderka, D., J. Le, and J. Vilecek. 1989. IL-6 inhibits lipopolysaccharide-induced tumor necrosis factor production in cultured human monocytes, U937 cells, and in mice. *J. Immunol.* **143**:3517-3523.
52. Snyman, E. and J. V. Van der Merwe. 1986. Endotoxin-polluted medium in a human in vitro fertilization program. *Fertility and Sterility* **46**:273-276.
53. Fishel, S., P. Jackson, J. Webster, and B. Faratian. 1988. Endotoxins in culture medium for human in vitro fertilization. *Fertilization and Sterility* **49**(1):108-111.
54. Orr, H. C., K. H. Sibinovic, P. G. Probst, and D. C. Littlejohn. 1975. Bacteriological activity in unfiltered calf sera collected for tissue culture use. *In Vitro* **11**:230-233.
55. Fumarola, D. and E. Jirillo. 1976. Contamination of con A preparations. *J. Immunol.* **116**:1197.
56. Fumarola, D. and E. Jirillo. 1979. Endotoxin contamination of some commercial preparations used in experimental research. *Biomedical Applications* **29**:379-385.
57. Del Llano, A. M. and J. A. Lavergne. 1991. Submitogenic doses of lipopolysaccharide alter the patterns of nuclear granularity and cell proliferation in human mixed lymphocyte cultures. *Transplantation Proceedings* **23**:1766-1770.
58. Wollman, Y., Y. Barak, and S. Levin. 1982. Contamination of human placental granulocyte macrophage colony stimulating factor by bacterial endotoxin: Its removal by a simple adsorption technique. *Biomedicine and Pharmacotherapy* **36**:323-325.

APRIL

April 27 - 29

LAL Methodology and Applications Seminar and Workshop

Holiday Inn-Crabtree
Raleigh, North Carolina

MAY

May 11 - 12

Pharmaceutical Water Monitoring and Testing

*LAL testing segment presented by Marilyn J. Gould, Ph.D.*PDA Training and Research Institute
Baltimore, Maryland
For more information, call (301) 986-0293.

May 18 - 20

ASM Annual Meeting

Georgia World Congress Center, Booth 916
Atlanta, Georgia

JUNE

June 15 - 17

LAL Methodology and Applications Seminar and Workshop

San Juan Marriott
San Juan, Puerto Rico

**For customer service, call
(800) LAL-TEST or
(508) 540-3444.**

**For technical service, call
(800) 848-3248 or
(508) 540-3444.**

**Please visit our website!
www.acciusa.com**

©1998 Associates of Cape Cod, Inc.
*All rights reserved.
Printed on recycled paper.*



Creating New Horizons in Endotoxin Testing

704 Main Street ■ Falmouth, MA 02540